




MaQS version 4.1.24 Release Notes

This document details the changes in this release of the MaQS product.

1. Selecting a Default Process for Sterilisation and Washing

In the Equipment Wizard under Cycle Combinations by selecting a process and clicking *Change*, the following pop-up screen will be displayed.




The process chosen can be selected as the 'Default' process for future Pack Template Processes. This means when creating a new pack template the information on the Process tab will be automatically updated with the default processes chosen.

These processes can still be changed if required by opening the Pack Template Process tab and selecting change.

2 New tab when printing labels

The print label screen has a new tab:




when printing a label by searching for the instrument. This separation from the previous version will remove any confusion when searching via an instrument.

3 Access to Help File.

The 'Help' menu on the MaQS client has been extended to allow access to the online help file and to the 'change history' document.

4 Change instrument quantity in Pack Templates

The Template Contents tab in the Pack templates screen now allows the quantity of an individual instrument template to be changed:



5 Retrieve instrument template details from the Pack Template screen

Double-clicking any instrument template in the 'Template contents' tab of the Pack Template screen will retrieve the details of that Instrument Template.

6 Indicator of no physical instruments

When an instrument template is added to the 'Template contents' tab of the Pack Template screen and the instrument template has no physical instruments, the entry will be shown in *italics*.

7 Installing software for the PDA/Handheld scanner

The 'Handheld/Scanner' menu now has an additional menu item 'Install PDA software' which will initiate the program for installing MaQS software on the scanner. The installation software is deployed to the CEInstaller subdirectory and is used to install MaQS programs on a scanner attached to the PC via Active Sync.

8 Validation of Operator IDs

The desktop scanning application will now ensure that an Operator ID has not been terminated.

9 Changing quantity in Preference Cards

When adding an item to a preference card and that item already exists, the existing item's quantity is incremented instead of creating a new entry in the preference card.

10 Consumable Packs

A Consumable Pack can now be deleted. However, if a subsequent allocation of that pack is made to a patient, it will result in an error.

The 'Remove All' option in Consumable Pack contents now works.

11 Cancel Patient Items

This screen now uses the current year and month as the default.

12 MDLookup Application


Columns in this application can now be sorted in ascending/descending order.

13 Default Quantity value to '1'

The handheld scanning program now sets the value to '1' when entering a quantity.

14 Inventory Management (Buy-Ins/Presterile)

Choose this option in System Settings by placing a tick in the box under Inventory Management if you do not want MaQS to manage inventory. Buy-ins/Pre-sterile items can still be issued to stores and patients without having to first receive items into the issuing store. Inventory reports will not be available if this option is selected.




15 Return Reasons and Checking Codes

Both of these options now allow a maximum of twenty Return Reasons and Checking Codes to be used. The scanning screens have been changed to allow the extra items to be shown/selected.

16 Where Is Patient information

The WhereIs screen has been modified to show the Items used on a patient:



17 Images and MaQS

Images of Packs and Instruments are stored in the MaQS database. Digital images are captured at various resolutions (640x480, 1024x768 etc) and the higher the resolution, the larger the size (in KB or MB) of the image, the more space it occupies in the database and the longer it takes to store into and retrieve from, the database.

For purposes of viewing on screen or for printing in a report, the higher resolutions have no added benefit. In some instances, users have reported being unable to print a report because the image it requires is too big.

When importing images, MaQS will now ensure that the larger images are resized to a resolution of 640x480.

For existing images in MaQS which may be too large, a utility is available which can resize them to the 640x480 resolution:



The utility shows the width and height (resolution) of the image and if it is higher than 640x480, the 'Resize!' button is enabled. Clicking the button will resize the image and update the database.

Contact MaQS Support if you would like a copy of the utility.